

ThirdEye^{SDI} this dental camera has convinced me!

Order Hereby we order the following products (price list 01/2025; prices in €uro)

Item-No	Item description	unit price	quantity	total
SDI_001*	ThirdEye^{SDI} dental camera set full-HD, 10x opt. zoom, autofocus, HD-SDI out, microphone, freeze, remote control <i>including:</i> miniature camera mount (with installation set) medical power supply (110/240volts AC → 12 volts DC) camera cable (HD-SDI/power/audio) white, 6m wireless remote control (handheld) Heliopan close-up lens #1 (mounted) + optional: close-up lens #2	3.480,00		
SDI_002 ¹	miniature camera mount (with M2 Allen screwdriver and Hylosil [®] -glue)	45,80		
SDI_003 ¹	medical power supply (110/240 volts AC → 12 volts DC)	105,00		
SDI_004 ¹	camera cable (HD-SDI/power/audio) white, 6m	165,00		
SDI_005	camera cable (HD-SDI/power/audio) white, extra meter	5,00		
SDI_006	mount adapter for Sirona C4/LED C8, XO-Flex and some other dental lights	55,00		
SDI_DVR	SDI recorder , recording format MP4 (5GB/h) SDI-in, microphone-in, HDMI-out, power supply, incl. 128 GB USB flash drive	295,00		
	¹ = equipment needed to use one single camera in two or more rooms			
SDI_011	Rental* of ThirdEye^{SDI} dental camera for testing: one month * Rental fee is fully credited, if purchasing camera !	345,00		
SDI_012	Rental* of ThirdEye^{SDI} dental camera for testing: two weeks * Rental fee is fully credited, if purchasing camera !	195,00		
	International camera rentals require the full sales price as a deposit. Deposit will be refunded, after camera was sent back w/o damage!	amount		
	* only to be paid for sales within European Community (if VAT-ID not provided)	+ 19 % VAT*		
		invoice total		

To ensure a correct delivery we need some more details:

Type of dental light room 1 room 2 (e.g. Sirona LED)

Length of camera cable room 1 m room 2 m

name	
address	
ZIP / city	
country	
telephone	facsimile

stamp

delivery in week 2025 (delivery time: 1-2 weeks) date signature

**All prices without costs for shipping, transport insurance and 19% German VAT* (*only within EU)
International sales against advance payment only. Please see our terms of business.**

Dr. Benno Raddatz Verlag Neue Medien Rheinstraße 29 76448 Durmersheim / Germany
Tel: +49 - 7245 - 937186 Mobile: +49 - 170 - 3812388 email: dr.raddatz@t-online.de